

Do you also want to administer the gentlest of care to a new-born baby?

The new-born patient needs all your professional skill and care – as well as lots of undisturbed rest. By choosing Unilect™ electrodes, you are helping to give the best possible care.

Unilect™ electrodes for the new-born

■ Are suited to their special needs and conditions

The special bio-adhesive gel ensures good adhesive properties that are suited to the humidity of an incubator – so electrodes don't need changing as frequently. The electrodes and cables are so soft that they cause no irritation, which means the baby is not unnecessarily disturbed - and neither is your work routine.

■ Do not affect adversely scanning and X-ray examinations

Radio-translucent electrodes such as Unilect™ are compatible with MRI scans and X-ray procedures.

■ Are easy to use






The electrodes are packed in bags of three, exactly the number you need for each patient. They are easy to unpack and apply and the gel does not stick to clothes and bed linen.





■ The electrodes have been specially developed to be gentle on the skin of newborn babies and to be easy to use.

Unilect™ Electrodes for Neonatal Department (Monitoring and 12-lead ECG)

Product	Ref. No.	Gel	Material	Connector & Cable length (mm)	Dimensions (mmxmm / mm2)	Packaging (liner/pouch/box/carton)
	40550	Bio-adhesive gel	Cloth	4 mm, 600 mm	22x27/594	3/3/300/3000
	40554	Bio-adhesive gel	Cloth	1.5 mm DIN, 600 mm	22x22/484	3/3/300/3000
	40555	Bio-adhesive gel	Cloth	1.5 mm DIN, 600 mm	22x27/594	3/3/300/3000
	40556	Bio-adhesive gel	Cloth	4 mm, 600 mm	22x22/484	3/3/300/3000
	40801	Bio-adhesive gel	Cloth	4 mm, 100 mm	19x26/494	5/10/300/3000

Product Data Sheets are available on request from your local Unomedical Sales Representative or at www.unilect-electrodes.com



Unilect™ Electrodes for Neonatal Department